Harvard Ionoglas Fill Extra

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7052000

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7092000

Conventional Glass Ionomer Cements

Harvard Ionoglas Fill Extra

Esthetic conventional glass ionomer cement for restorations, self cure

- Esthetic natural translucency
- High surface hardness, reliable mechanical properties
- Low solubility
- Very good marginal seal
- Excellent bond to dentin and enamel without etching
- High fluoride release
- Radiopaque
- Also available in OptiCaps®

High esthetics!

Harvard Ionoglas Cem Extra

Esthetic conventional glass ionomer cement for luting of crowns, bridges, inlays and onlays

- Perfect consistency and low film thickness
- Very good marginal fit and seal
- Excellent bond to dentin and enamel without etching
- Reduced solubility
- High fluoride release
- Also available in OptiCaps®

Harvard Ionoglas Fill Extra

Properties and Advantages!

OptiCaps® and Handmix

- *Reliable and convenient*
- *Easy handling*

Nice translucency and very easy to polish

- *Esthetic final results*

Excellent bond to dentin and enamel without etching

- *Very good marginal fit and seal*

High surface hardness


Low solubility

No polymerization shrinkage

High fluoride release

Radiopaque

High fluoride release

Cumulative fluoride release*

Application


Prepared cavity


Finished restoration

Good physical properties

Flexural strength*

Compressive strength*

■ Harvard Ionoglas Fill Extra ■ GC Fuji IX GP (GC) ■ Ketac Fil Plus (3M Espe)

Initial protection:

Harvard Ionocoat LC

Prevents the glass ionomer restoration from early dissolution during setting.

Harvard Ionoglas Cem Extra

Properties and Advantages!

OptiCaps® and Handmix

- *Reliable and convenient*
- *Easy handling*

Excellent bond to dentin and enamel without etching

- *Very good marginal fit and seal*

Low shrinkage
Low solubility

High fluoride release

Radiopaque


Complementary to Harvard Ionoglas Fill Extra and Harvard Ionoglas Cem Extra:

Harvard Ionoconditioner


Removes the smear layer from the prepared tooth surface.

High fluoride release

Cumulative fluoride release*

Good physical properties

Flexural strength*

* Internal data by Harvard Dental International

Technical Data

Harvard Ionoglas Fill Extra

Esthetic conventional glass ionomer cement for restorations, self cure

Compressive strength	184 MPa
Flexural strength	34 MPa
Working time	1:30 min
Setting time	3:30 min
OptiCaps® mixing time	10 sec
Working time	90 sec (1:30 min)*

* from the start of mixing

Harvard Ionoglas Cem Extra

Esthetic conventional glass ionomer cement for luting of crowns, bridges, inlays and onlays

Compressive strength	92 MPa
Flexural strength	16 MPa
Working time	1:30 min
Setting time	5:00 min
OptiCaps® mixing time	10 sec
Working time	90 sec (1:30 min)*

Article

Order no.

Harvard Ionoglas Fill Extra

15 g powder / 8 ml liquid, dosage spoon, mixing pad

Shade A2	7052112
Shade A3	7052113
Shade A3.5	7052135

50 OptiCaps® ea. 0.5 g

Shade A2	7052252
Shade A3	7052253
Shade A3.5	7052254

Harvard Ionoglas Cem Extra

15 g powder, shade Universal / 10 ml liquid, dosage spoon, mixing pad

50 OptiCaps® ea. 0.4 g, shade Universal	7042115
	7042250

Harvard Ionocoat LC 5 ml bottle	7052000
---------------------------------	---------

Harvard Ionoconditioner 5 ml bottle	7053000
-------------------------------------	---------

Harvard Applier OptiCaps®	7092000
---------------------------	---------

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V02.1 11/2016
10331

Harvard Ionoglas Cem Extra

Classical Glass Ionomer Luting Cement

Harvard Ionoglas Cem Extra is a new glassionomer luting cement with improved physical properties.

Beside its high fluoride release and the excellent biocompatibility as glassionomer cement Harvard Ionoglas Cem Extra has also good chemical bonding to dentine and enamel and to the restoration and a tight seal at the dentinal margins. Because of its radiopacity it ensures easy postoperative diagnosis.

Indications

- cementing of crowns and bridges (all types: metal, ceramic fused to metal and ceramic)
- cementation of metal inlays and onlays
- liner under composite fillings

Contraindications / side-effects

Do not use Harvard Ionoglas Cem Extra with patients who show an allergy to the material. If allergic reactions occur, further treatment with Harvard Ionoglas Cem Extra must be discontinued and the patient must be advised to consult a physician.

Do not allow the liquid or cement mixture to contact the oral tissues or skin. In case of contact, remove the material with absorbent cotton soaked in alcohol and rinse with water.

Avoid eye contact of the liquid or cement mixture. In case of contact, immediately flush with water and seek medical advise.

Interaction with other materials

Do not mix the powder or liquid of Harvard Ionoglas Cem Extra with any other glass-ionomer product.

Application

1. Preparation

Prepare the tooth to be treated in the accustomed manner.

Clean the prepared tooth with pumice and water. Rinse thoroughly and dry, but **do not** desiccate.

By using a cotton pellet apply a dentine conditioner (e.g. Harvard Ionoconditioner) for **20 sec** to remove the smear layer. Rinse the dentine conditioner with water and dry it in an airstream, but **do not** desiccate.

Pulp capping with **Harvard Ionoglas Cem Extra** is **contraindicated**. To deep areas of possible pulpal exposure apply a small amount of a calcium hydroxide liner (e.g. Harvard CalciumHydroxide).

2. Dosing of powder and liquid

The powder/liquid ratio to achieve a suitable consistency is 1.8 / 1.0. Mix 1 level spoon (blue) of powder and 2 drops of liquid.

Shake the bottle to loosen the powder. Overfill the spoon with the powder, level the powder for using the collar at the top of the bottle and carry it onto the mixing pad provided. Avoid compressing powder into the spoon with the inside wall of the bottle.

Turn the liquid bottle vertically with the tip about 5cm above the mixing pad. Steady your hand and squeeze the bottle gently to dispense one drop. If any bubbles are present, lightly tap the bottle with the fingers holding it. **Discount drops** that are obviously not full-sized.

3. Mixing of powder and liquid

Use a small spatula to rapidly mix all the cement powder into the liquid. Usually one scoop powder/ two drops liquid should provide for sufficient amount of mixed cement. The mixed cement should be thixotropic and have a smooth consistency and glossy appearance. Total mixing time is **30 sec**.

After use, tightly close both liquid and powder bottles to prevent exposure to moisture.

4. Cementation Technique

Mix the required amount of **Harvard Ionoglas Cem Extra**. Coat the prepared bonding surface of the restoration with **1 mm** of cement (**Do not overfill!**) and seat immediately **within the working time**. Seat the restoration by using steady, even pressure.

Remove excess cement at the first setting stage.

Maintain isolation until the set of the cement is verified.

Working time from start of mix at 23°C (74°F): **1:30 min**
Net setting time at 37°C (99°F): **5:00 min**

Notes

Higher temperatures will shorten the working time, lower temperatures will prolong the working time.

An overextended working time reduces adhesion to enamel and dentine.

Storage information

Store **Harvard Ionoglas Cem Extra** in a cool place at temperatures between 4 and 25° C or 39 and 77°F!

Do not use after expiry date.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and Harvard Dental Internationals sole obligation shall be repair or replacement of the Harvard Dental International product.

Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Only for dental use!

Store product out of reach of children!

Harvard Ionoglas Cem Extra conforms to the DIN EN ISO 9917-1.

The instructions for use has to be kept for the duration of the application.

Order-No. Article

7042115 15 g powder, shade Universal / 10 ml liquid,
dosage spoon, mixing pad

 Manufacturer	Harvard Dental International GmbH Margaretenstr. 2 – 4 15366 Hoppegarten, Germany	 0482 Made in Germany	Tel: +49 (0) 30/99 28 978-0 Fax: +49 (0) 30/99 28 978-19 info@harvard-dental-international.de www.harvard-dental.de	V02.1 11/2016 10331
---	--	--	--	------------------------