Harvard PremiumSil

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Harvard PremiumSil

Super *hydrophilic* impression materials

- Super hydrophilic for perfect wetting and perfect accuracy in all clinical situations
- Long working time and short setting time
- Contrasting colors – for optimal readability
- Pleasant berry smell
- Automatic mixing also for putty

Because everything begins with a good impression!

Precise dental restorations begin with a good impression. Dental impression materials are the key for clinical success.

Harvard presents a new series of A-silicone (vinylpolysiloxane) impression materials: **Harvard PremiumSil** products.

The Harvard A-Silicones combine the best properties:

- **excellent flow properties** to minimize voids
- **high hydrophilicity** for excellent reproduction of details
- **long working time** and **short setting time** comfortable for dentist and patient
- very low **shrinkage**
- excellent tear strength and elasticity for **distortion-free results**
- optimal **kneadable consistency** for the putty products
- smart **automatic mixing** also for putty
- **range of products** for all clinical needs
- pleasant smell

Technical data

All **Harvard PremiumSil** products are offered in two setting times. Therefore each practitioner will find the perfect material for his specific clinical need. Due to snap set technology the setting time is minimized, which means less waiting time for patient and dentist.

The snap set technology of **Harvard PremiumSil** also leads to excellent elastic recovery. Elastic recovery is the ability of the impression to return to its original shape after mechanical stress, for instance after removal of undercuts. Thanks to the perfect snap set technology, **Harvard PremiumSil** is completely cured and has a very high degree of elasticity when it is removed from the mouth.

All **Harvard PremiumSil** products have an elastic recovery of > 99,5 %.

With **Harvard PremiumSil**, high quality in all material parameters is the first priority.

Harvard PremiumSil products*	Working time (min:sec)	Time in mouth (min:sec)	Setting time (min:sec)	Shore A hardness after 24 h
Light	01:45	02:30	03:30	50
Medium				56
Mono				62
Heavy	01:15	01:30	02:30	65
Putty Soft				69
Putty				71
Putty Soft <i>MaxiMix</i>	00:50	01:45	02:45	68
Light Fast				47
Medium Fast				57
Mono Fast	01:15	01:30	02:30	63
Heavy Fast				64
Putty Soft Fast				68
Putty Fast	01:15	01:30	02:30	70


* Internal data by Harvard Dental International

Snap set properties

A comfortable working time adjusted to the needs of the practitioner is very important.


After placement the material has to set. The faster this setting, the more convenient it is for the patient. A long working time with a short setting time is called: „snap set“.

Comparison of **Harvard PremiumSil Light** and **Harvard PremiumSil Light Fast**

2
setting
times

Comparison with products


■ Harvard PremiumSil Light ■ Aquasil Ultra XLV (Dentsply)
■ Flexitime Correct Flow (Heraeus Kulzer) ■ Honigum Light (DMG)

* Internal data by Harvard Dental International

Super hydrophilic impressions

Harvard PremiumSil displaces moisture and forms an intimate contact with the tooth and tissue surface. Reliable for perfectly fitting restorations.

To measure this property, a drop of water is placed on the impression material and the spreading on the surface is observed.

	immediately after placement	after 5 sec.	after 10 sec.	after 30 sec.	after 60 sec.
Harvard PremiumSil Light					
Elite HD+ (Zhermack)					
Aquasil Ultra LV (Dentsply)					

Harvard PremiumSil super hydrophilic properties are optimal to capture the highest degree of surface details by excellent wetting.

Overview according to technique

Harvard PremiumSil product portfolio

2
setting
times

Harvard PremiumSil Light
Vinylpolysiloxane impression material,
super hydrophilic, light bodied, type 3

Harvard PremiumSil Medium
Vinylpolysiloxane impression material,
super hydrophilic, medium bodied, type 3

Harvard PremiumSil Mono
Vinylpolysiloxane impression material,
super hydrophilic, monophasic, type 2

Harvard PremiumSil Heavy
Vinylpolysiloxane impression material,
super hydrophilic, heavy bodied, type 2


Harvard PremiumSil Putty Soft *MaxiMix*
Vinylpolysiloxane impression material,
Putty consistency, type 0

Harvard PremiumSil Putty Soft
Vinylpolysiloxane impression material,
Putty Soft, kneadable, type 0

Harvard PremiumSil Putty
Vinylpolysiloxane impression material,
Putty, kneadable, type 0

Harvard TraySive
Tray adhesive

Very berry! Snappy set!

The recommended combinations of all **Harvard PremiumSil** products according to the used techniques:

Techniques	Harvard PremiumSil viscosity combinations			
2 steps 2 phases (putty-wash technique)	PremiumSil Putty + PremiumSil Light	or	PremiumSil Putty + PremiumSil Medium	
1 step 2 phases (sandwich technique)	PremiumSil Putty Soft + PremiumSil Medium	or	PremiumSil Heavy + PremiumSil Medium	
1 step 1 phase	PremiumSil Mono			

Harvard PremiumSil Putty Soft *MaxiMix*:
Automatic mixing also for putty - easy and moldable!

Harvard TraySive
Tray adhesive

Article	Color	Order no.
2 automix cartridges ea. 50 ml, 12 mixing tips 1:1: Harvard PremiumSil Light Harvard PremiumSil Light Fast Harvard PremiumSil Medium Harvard PremiumSil Medium Fast Harvard PremiumSil Mono Harvard PremiumSil Mono Fast Harvard PremiumSil Heavy Harvard PremiumSil Heavy Fast	orange orange bright green bright green bright blue bright blue bright purple bright purple	7083701 7083702 7083711 7083712 7083721 7083722 7083731 7083732
380 ml maximix cartridge, 10 dynamic mixers, 1 bayonet ring: Harvard PremiumSil Putty Soft MaxiMix, type 0 Harvard PremiumSil Heavy Fast MaxiMix, type 2	purple bright purple	7083740 7083733
2 jars ea. 300 g, 2 spoons: Harvard PremiumSil Putty Fast, type 0	dark purple	7083763
2 jars ea. 600 g, 2 spoons: Harvard PremiumSil Putty Harvard PremiumSil Putty Fast Harvard PremiumSil Putty Soft Harvard PremiumSil Putty Soft Fast	dark purple dark purple purple purple	7083751 7083752 7083741 7083742
Harvard Auto 1:1 O-Green, refill bag with 50 mixing tips for Harvard PremiumSil Heavy / Heavy Fast		7096000
Harvard Auto 1:1 O-Pink, refill bag with 50 mixing tips for Harvard PremiumSil Mono / Mono Fast		7091200
Harvard Auto 1:1 O-Yellow, refill bag with 50 mixing tips for Harvard PremiumSil Light / Light Fast; Harvard PremiumSil Medium / Medium Fast		7091300
Harvard IntraTips yellow Refill bag with 50 IntraTips for Harvard PremiumSil Light Harvard PremiumSil Light Fast		7083619
Harvard Maxi 5:1 Dynamic, refill bag with 50 dynamic mixers and 2 bayonet rings for Harvard PremiumSil Putty Soft MaxiMix; Harvard PremiumSil Heavy Fast MaxiMix		7091400
Harvard Dispenser Automix 1:1		7095100
Harvard TraySive 10 ml bottle		7083753

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Instructions for Use EN

Vinylpolysiloxane impression material, super hydrophilic

Harvard PremiumSil Medium is very versatile. It can be used:

- as a hydrophilic syringeable material for the (one step) two phase impression technique in combination with Harvard PremiumSil Putty Soft or Harvard PremiumSil Heavy
- as a correction material in the (two step) putty/wash technique in combination with Harvard PremiumSil Putty
- as a syringable one phase impression material for crowns, bridges, inlays, partial and complete dentures. It also serves as an impression material for relining, functional impressions, impressions for denture repairs and for impressions using the copper band technique.

Harvard PremiumSil Medium is a light to medium viscosity hydrophilic precision impression material based on addition curing vinyl silicones delivered in automix-cartridges. It can be extruded directly onto the teeth or the tray.

Harvard PremiumSil Medium shows strong hydrophilic properties. Therefore it improves the wettability and yields an exact reproduction of subgingival areas and other moist surfaces.

Harvard PremiumSil Medium is thixotropic. That means it stays in the position in the tray and on the tooth, but under pressure it flows.

Harvard PremiumSil Medium meets the requirements of ISO 4823, type 3 (light-bodied consistency).

Contraindications/side effects

None known

Application

1. Preparation of the automix-cartridge

Harvard Dental International GmbH only recommends for Harvard PremiumSil Medium the use of mixing tips type Harvard Auto 1:1 O-Yellow.

Mount the cartridge into the application gun. Remove and dispose the closure-cap. For exact flow control extrude slightly material until uniform flow from both orifices is achieved. Install the mixing tip onto the cartridge and turn it 90° clockwise until it locks in position, and extrude needed material.

For direct extruding onto teeth use Intra Tips yellow.

Note:

Store used cartridge with fixed used mixing tip. **Do not re-use the cartridge cap.** Closing the cartridge by turning on the cap between use could cause cross contamination of components and thus premature setting.

Leave the used mixing tip on the cartridge. It serves as a cap.

2. Impression taking

Load the auto mixed material onto the tray or directly onto the teeth. Insert the tray not later than **1:10 min** (working time at 23°C or 74°F) after extruding Harvard PremiumSil Medium into the mouth. Maintain position for at least **1:30 minutes**.

3. Casting

The impression should be stored under normal room conditions (max. 25°C or 77°F). The impression can be poured 30 minutes after removal from the mouth. No other time limitation does exist.

4. Electroplating

All impressions can be copper- or silverplated.

5. Disinfection

Impressions can be disinfected with standard disinfection solutions.

Notes:

Do not contact the product with latex gloves, they may inhibit the setting. Do not touch teeth with latex gloves.

Do not use with standard (condensation reaction type) silicone impression materials.

Storage

Do not store above 25 °C (77 °F). Do not use after expiry date.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only!

Store product out of reach of children!

Keep the instructions for use for the duration of the application.

Order-No.	Article
7083712	Harvard PremiumSil Medium Fast, 2 automix cartridges ea. 50 ml, 12 mixing tips 1:1, 4 Intra Tips
7091300	Harvard Auto 1:1 O-Yellow, refill bag with 50 mixing tips
7083619	Harvard Intra Tips yellow refill bag with 50 intra tips.
7095100	Harvard Dispenser Automix 1:1

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