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7051000

7092000

Resin modified glass ionomer cements

Harvard Ionoresin Fill

Light cure esthetic glass ionomer cement
for restorations

- Ultra fine fillers
- Nice polishability
- Very low shrinkage
- Very good marginal fit and seal
- Three cure mode:
self cure + light cure + classical cement reaction
- Also available in OptiCaps[®]

Harvard Ionoresin Cem

Light cure resin modified glass ionomer cement
for luting and lining

- Ultra fine fillers
- Low film thickness
- Very low shrinkage
- Very good marginal fit and seal
- Three cure modes:
self cure + light cure + classical cement reaction
- Also available in OptiCaps[®]

Harvard Ionoresin Fill

Properties and Advantages!

OptiCaps® and Handmix

- *Reliable and convenient*
- *Easy handling*

Ultra fine fillers

very easy to polish

- *Esthetic final results*

Three cure mode:

self cure + light cure + classical cement reaction

- *Cures also in dark areas*

Bonds to dentin

- *Very good marginal fit and seal*

Higher mechanical strength than conventional GICs

Virtually insoluble


Higher surface hardness than conventional GICs

Very low shrinkage

Strong and long-lasting esthetic restorations

Radiopaque

Good physical properties

Application

Prepared cavity

Finished restoration

Recommendation:

Harvard Ionoresin Prime LC

Optional primer for enamel and dentin for best results

Harvard Ionoresin Cem

Properties and Advantages!

OptiCaps® and Handmix

- *Reliable and convenient*
- *Easy handling*

Ultra fine fillers

- *Low film thickness*

Three cure mode:

self cure + light cure + classical cement reaction

- *Very good marginal fit and seal*

Excellent bond to dentin

Higher mechanical strength than conventional GICs

Virtually insoluble

Very low shrinkage

Higher surface hardness than conventional GICs


Strong and long-lasting restorations

Radiopaque


save time with OptiCaps®

Good physical properties

Flexural strength*

Compressive strength*

* Internal data by Harvard Dental International

Technical Data

Harvard Ionoresin Fill

Resin modified glass ionomer cement for restorations, light cure

Compressive strength - self-curing	220 MPa
Compressive strength - light-curing	250 MPa
Flexural strength	80 MPa
Barcol hardness	72

Working time	2:00 min
Setting time (self-curing)	4:30 min
Setting time (Light-curing)	20 sec

Harvard Ionoresin Cem

Resin modified glass ionomer cement for luting and lining, light cure

Compressive strength - self-curing	190 MPa
Compressive strength - light-curing	214 MPa
Flexural strength	70 MPa

Working time	2:30 min
Setting time (self-curing)	4:30 min
Setting time (Light-curing)	20 sec

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Users Instruction EN

Harvard Ionoresin Cem**Harvard Ionoresin Cem OptiCaps®**

Resin modified glass ionomer cement for luting and lining, light-curing

Harvard Ionoresin Cem is a light-curing, resin modified glass ionomer luting and lining cement.

Good bonding characteristics between dentine, enamel and restoration assure a good marginal fit. **Harvard Ionoresin Cem** is biocompatible, X-ray opaque and releases fluoride ions.

Harvard Ionoresin Cem OptiCaps® can quickly and simply be activated and extruded with the Capsule applicer. The mixing of the capsule is done in 10 s by means of a common high frequency mixer at 4,300 oscillations/min. The cement can be applied directly from the capsule.

Indications

- Cementing of crowns, bridges, inlays and onlays (all types: metal, resin, ceramic fused to metal and ceramic)
- Liner under composite fillings

Composition

Barium aluminium silicate glass, polyacrylic acids, methacrylates

Contraindications / side-effects

Harvard Ionoresin Cem may not be used with patients with a hypersensitivity to any of the components. If allergic reactions occur, further treatment with **Harvard Ionoresin Cem** must be discontinued and the patient must be advised to consult a physician.

The contact of **Harvard Ionoresin Cem** liquid or mixed cement with the skin or oral tissue should be avoided. In case of accidental contact, remove the material with alcohol soaked cotton pads. Rinse the moistened places with plenty of water.

In case of contact with the eye, rinse with plenty of water with the eyelid opened and consult an ophthalmologist.

Interaction with other materials

The mixing of powder and liquid with components of other glass ionomer cements is not allowed.

Avoid direct contact with products containing eugenol since eugenol can impair the curing of **Harvard Ionoresin Cem**.

Application**1. Preparation**

Prepare the tooth to be treated in the accustomed manner.

Using a micro brush, apply **Harvard Ionoresin Prime LC** onto the moist dentine and enamel surfaces for 30 s. Note that the material should be distributed homogeneously over the entire area. By gently blowing with compressed air (10 s) the volatile components are removed and the primer uniformly dispersed. Subsequently, light-cure the layer with a suitable dental curing light (10 s).

To areas of possible pulpal exposure apply a small quantity of calcium hydroxide liner, as the capping of the pulp is contraindicated.

2. Dosing of powder and liquid (only for hand mix)

The powder to liquid ratio amounts to 1.7 / 1.0.

Mix 1 level spoon (blue) of powder and 3 drops of liquid.

For accurate dispensing of **Harvard Ionoresin Cem** powder shake the bottle to loosen the powder. Overfill the spoon with the powder, level the powder using the scraper at the top of the bottle and carry it onto the mixing pad. Avoid compressing powder into the spoon at the inside wall of the bottle.

The bottle with the liquid is held vertically about 2 cm over the mixing pad. Using slight pressure, the liquid is dropped next to the cement. For mixing, use only drops without bubbles.

3. Mixing**3.1 Mixing of powder and liquid (hand mix)**

With a plastic spatula, the entire powder quantity is mixed within 30 s into the liquid. After mixing, the cement must display a homogeneous consistency and thixotropic behaviour.

After use, firmly close both bottles in order to prevent exposure to moisture!

3.2 Mixing Harvard Ionoresin Cem OptiCaps®

The activation and mixing of the OptiCaps® is done corresponding to the directions of the "Instructions for the activation and mixing of Harvard OptiCaps®".

The mixing time is 10 s.

Avoid lag-times between the processes of activation, mixing and application, as the material is in the process of setting and time-lags may impair or prevent application of the material.

Hence, all the cement required should be dispensed from the Harvard OptiCaps® within 30 s from the start of mix.

4. Cementing

The restoration is prepared in accordance with manufacturer instructions.

4.1 Hand mix

Mix the required quantity of **Harvard Ionoresin Cem** and apply to the surface areas of the restoration in a thin layer (max. 1 mm in height). Insert the restoration immediately.

Working time from start of mix at 23°C or 74°F:

2:30 min

Net setting time without light-curing at 37°C or 99°F:

4:30 min**4.2 Cementing when using Harvard OptiCaps®**

Apply a layer (max. 1 mm high) of **Harvard Ionoresin Cem** to the surface areas of the restoration and insert the restoration immediately.

Working time from start of mix at 23°C or 74°F:

2:00 min

Net setting time without light-curing at 37°C or 99°F:

4:30 min**4.3 References for hand mix and OptiCaps®**

Higher temperatures shorten, lower temperatures lengthen the working time of the cement.

Processing the restoration for longer than the working time of the cement diminishes the adhesion of the cement to enamel and dentine.

Large cement excesses can be removed during the setting phase.

Hold the tooth isolated until the cement is set.

In order to reach an optimal setting, the cement layer should be light-cured for 20 s.

Storage information

Store **Harvard Ionoresin Cem** in a cool and dark place at temperatures between 4 and 25°C or 39 and 77 °F!

Do not use the product beyond the end of the expiration date!

Only for dental use!

Store product out of reach of children!

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

Harvard OptiCaps[®] before activation
The capsule is shaken to loosen the powder easily

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To activate the Harvard OptiCaps[®] totally depress the plunger on a hard and plain surface.

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
Insert the Harvard OptiCaps[®] into the capsule applicator and click once to standardize. Note: The plunger must be at the same level as the bottom of the capsule

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
Insert the Harvard OptiCaps[®] into a capsule mixer. Close lid and mix immediately for **10 s** at 4300 oscillations / min.

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After the end of mixing remove the pin from the nozzle. Otherwise capsule can burst.

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Insert the Harvard OptiCaps[®] into the capsule applicator. Squeeze lever 2 times (2 clicks) to prime the Harvard OptiCaps. Extrude the mixed material directly into the preparation. Unlock the gun (push button A) and remove the Harvard OptiCaps.