HARVARD®

# Harvard

## MTA Universal Maišomas rankomis

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# Harvard

## MTA Universal

Endodontic **biocompatible** repair cement

## HandMix

HandMix:  
flexible &  
economical

## Properties

- Highly biocompatible material
- Very good bonding to dentin
- No solubility in tissue fluid after setting
- High pH-value
- HandMix: flexible and economical
- Radiopaque
- Next clinical step after 5 minutes

## Advantages

- *Repairing of tooth structure*
- *Formation of tertiary dentin over the pulp ("dentin bridge")*
- *Formation of a protective layer against bacteria due to the formation of calcium hydroxide and apatite*
- *Setting unaffected by humidity*
- *Can easily be mixed to a nice creamy consistency*

# One biocompatible material: several indications!

## Overview of the clinical indications of Harvard MTA Universal HandMix

### Pulp Capping

- Pulp protection:  
direct and indirect pulp capping

### Filling of root canal perforation


- Excellent seal
- Formation of a protective layer preventing bacterial penetration
- Biocompatible

### Root end filling (orthograde)

- Sealing the apex in one step
- Compatible with guttapercha and all root canal sealers

### Root end filling (retrograde)

- Setting is not affected by humidity

## Technical Data

Endodontic **biocompatible** repair cement based on mineral oxides  
(mineral trioxide aggregate) as HandMix.

Mixing ratio Powder : Liquid (by weight)	2 : 1	2.6 : 1	3 : 1
Working time (23 °C / 73 °F)	3 min	2 min	1 min
Next clinical step after	5 min	5 min	5 min

## Article

## Order no.

### Harvard MTA Universal HandMix

1 g bottle powder, 3 ml bottle liquid,  
Dosing spoon, mixing pad, spatula

7081507

Also available:

### Harvard MTA Universal OptiCaps®

2 capsules ea. 0.3 g  
ea. packed in an aluminium pouch

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### 2. Application

#### 2.1. Repair of root perforations

Place rubber dam and clean the root canal system using intracanal instruments and irrigate with NaOCl. Dry the root canal with paper points and isolate the perforation. Fill the root canal apical of the perforation completely with a suitable root canal filling material. Mix Harvard MTA as described under point 1.

Apply Harvard MTA with suitable instruments into the perforation site and condense it.

Check the position of Harvard MTA in the root canal by an X-ray. If an adequate barrier has not been created, rinse the MTA cement out of the canal and repeat the procedure.

Remove excess moisture with a cotton pellet or a paper point.

Place a moist cotton pellet in the access to the root canal and apply a temporary filling material.

Alternatively seal the root canal with a suitable root canal filling material and seal the cavity with a tight filling. Both options can be done at the earliest **5 minutes** after placement of the Harvard MTA.

Harvard MTA repair material remains as a permanent part of the root canal filling.

### Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

### Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

### For dental use only!

#### Store product out of reach of children!

Keep the instructions for use for the duration of the application.

#### 2.2. Root end filling (retrograde)

Create a surgical access to the root-end and resect the root. Prepare an apical cavity to a depth of 3-5 mm.

Isolate the area and dry the root-end cavity with paper points. Achieve hemostasis with suitable methods.

Mix Harvard MTA as described under point 1.

Apply Harvard MTA with suitable instruments and condense it using a small plunger.

Remove excess cement and clean the surface of the root with a moist piece of gauze.

Confirm placement of Harvard MTA repair material with an X-ray. Harvard MTA remains as a permanent part of the root canal filling.

#### 2.3. Pulp Capping

Place rubber dam and prepare the cavity. Rinse the cavity and exposed pulpal areas with a suitable disinfectant.

Mix Harvard MTA as described under point 1.

Apply a small amount of Harvard MTA with a suitable instrument over the exposed pulp and remove excess moisture with a cotton pellet.

At the earliest **5 minutes** after application of Harvard MTA place a small amount of a flowable light cure liner (e.g. Harvard IonLine) and light cure.

Etch the remaining cavity walls according to the total-etch-technique with Harvard Etch and apply a suitable bonding agent (e.g. Harvard Bond TE Bond) according to the corresponding instructions.

Place a light cure composite (e.g. Harvard PremiumFill) according to the instructions and light cure.

Check pulp vitality and status regularly.

#### 2.4. Apexitifikation / Root-end filling (orthograde)

Place rubber dam and clean the root canal system using intracanal instruments and irrigate with NaOCl. Dry the root canal with paper points.

For disinfection place calcium hydroxide paste in the root canal and seal the access opening with a temporary filling material.

Remove the calcium hydroxide paste after one week from the root canal system.

Mix Harvard MTA as described under point 1.

With a suitable instrument apply a small amount of Harvard MTA into the apical region and condense it. Create a 3 - 5 mm barrier of MTA cement.

Check the position of Harvard MTA by an X-ray.

If an adequate barrier has not been created, rinse the MTA cement out of the canal and repeat the procedure.

At the earliest **5 minutes** after application of the Harvard MTA place a cotton pellet in the access to the root canal and apply a temporary filling material.

Alternatively seal the root canal with a suitable root canal filling material and seal the cavity with a tight filling.

Harvard MTA remains as a permanent part of the root canal filling.

#### Additional remarks

• In the first hour after application handle the placed MTA cement carefully.

• Intraloral application of Harvard MTA must be done immediately after mixing to prevent dehydration during setting.

• Harvard MTA can cause discoloration.

#### Storage

Store Harvard MTA at a dry place at 10 - 25 °C (50 °F - 77 °F). **Do not store below 10°C (50 °F)!** Do not use after expiry date.

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### Order-No.

7081507 Harvard MTA (Handmix), 1g powder, 3ml liquid, dosage spoon

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