Harvard Ionoglas Cem LC

Light cure resin modified glass ionomer cement for luting crowns and bridges

Properties

- Perfect consistency for cementation
- Fine fillers
- Light cure
- Adhesion to enamel and dentin
- Radiopaque
- Not sensitive to moisture

Advantages

- Durable restorations
- Easy to mix
- Convenient to apply
- Good fit and seal
- Esthetic restorations

Application

Harvard Ionoglas Cem LC is the perfect choice for the cementation of all kinds of restorations. Inlays, onlays, crowns and bridges made of metal, zirconia or ceramic can be luted with long-term clinical success.


Harvard Ionoglas Cem LC has improved mechanical properties compared to classical GIC. It is not sensitive to moisture. The application is very easy.

Harvard Ionoglas Cem LC results in esthetic and long-lasting luting and lining applications.

For optimal results use Ionoresin Prime LC.

Physical Properties

Flexural strength*

Compressive strength*

- Harvard Ionoglas Cem LC (self curing)
- ▨ Harvard Ionoglas Cem LC (light curing)
- GC Fuji Plus EWT (self curing)

* Internal data by Harvard Dental International

Technical Data

Light cure resin modified glass ionomer cement for luting and lining

Compressive strength - self curing	150 MPa
Compressive strength - light curing	160 MPa
Flexural strength	60 MPa
Working time	1:30 min
Setting time (self curing)	5:30 min
Light curing (optional)	20 sec

Article

Harvard Ionoglas Cem LC

15 g powder (shade Universal), 10 ml liquid, dosage spoon, mixing pad

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Harvard Ionoglas Cem LC

Resin modified glass ionomer cement for luting and lining, light-curing

Harvard Ionoglas Cem LC is a light-curing, resin modified glass ionomer cement for luting and lining.

As a glass ionomer cement it is simple to use. The curing is initiated by common visible light. Harvard Ionoglas Cem LC exhibits strong chemical adhesion bond to teeth and good biocompatibility. This glass ionomer cement is characterised by good mechanical properties.

Indications

- Cementing of crowns, bridges, inlays and onlays (all types: metal, resin, ceramic fused to metal and ceramic)
- Liner under composite fillings

Contraindications / side-effects

In singular cases, Harvard Ionoglas Cem LC may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used.

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of calciumhydroxide material.

Interaction with other materials

Avoid direct contact with products containing eugenol since eugenol impairs the setting of Harvard Ionoglas Cem LC.

Do not mix the powder or liquid of Harvard Ionoglas Cem LC with any other glass ionomer product.

Application

1. Preparation

Prepare the tooth to be restored in the usual manner.

For pulp protection areas close to the pulp should be covered with a thin layer of calciumhydroxide material (e.g. Harvard CalciumHydroxide) or MTA cement (e.g. Harvard MTA Universal).

Apply Harvard Ionoresin Prime LC with a suitable brush onto the **moist** enamel and dentine surfaces for **30 seconds** with agitation. The material should build a homogeneous layer. Air thin gently (**10 seconds**) to remove the volatile components and to disperse the adhesive. Then light cure with a suitable dental light cure unit (wavelength range 400–500 nm, light intensity min. 1000 mW/cm²) for **10 seconds**.

2. Dosing of powder and liquid

The powder to liquid ratio is 1.8 / 1.0. This can be obtained by mixing 1 level (orange) scoop of powder and 2 drops of liquid.

Note:

The spoon dosage is only an approximate reproduction of the nominal mixing ratio.

For accurate dispensing of Harvard Ionoglas Cem LC powder shake the bottle to loosen the powder. Overfill the spoon with the powder and level the powder using the scraper at the top of the bottle. Place the powder on the mixing pad.

For dispensing of Harvard Ionoglas Cem LC liquid turn the bottle vertically with the tip about 2 cm above the mixing pad. Squeeze the bottle gently to dispense one drop. If any bubbles are present, lightly tap the bottle with the fingers holding it. Discard drops that contain bubbles and are obviously not full-sized or oversized.

3. Mixing of powder and liquid

Use a small spatula to rapidly mix all the cement powder into the liquid. Usually one scoop powder/ 2 drops liquid should provide for sufficient amount of mixed cement. The mixed cement should be thixotropic and have a smooth consistency and glossy appearance. Total mixing time is **30 seconds**.

After use, tightly close both liquid and powder bottles to prevent exposure to moisture.

4. Cementing

Prepare the restoration according to manufacturer instructions.

Mix the required amount of Harvard Ionoglas Cem LC and apply a thin layer of max.1 mm to the prepared bonding surface of the restoration and seat immediately.

Working time from start of mix at 23°C (74°F): **1:30 min**

Net setting time without light curing at 37°C (99°F): **5:30 min**

Notes

Higher temperatures will shorten the working time, lower temperatures will prolong the working time.

An overextended working time reduces adhesion to enamel and dentine.

Remove excess cement at the first setting stage.

Maintain isolation until the set of the cement is verified.

Setting can be optimized with a **20 seconds** light cure.

Warnings

- Unpolymerized material may have an irritant effect and can lead to sensitization against methacrylates.
- Avoid contact with skin, mucous membrane and eyes
- If the material comes into contact with skin, In case of contact, remove the material with absorbent cotton soaked in alcohol and rinse with water. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.

Storage information

Storage temperature 4-25°C (39-77°F). Protect from direct sunlight. Do not use after expiry date.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only!

Store product out of reach of children!

Harvard Ionoglas Cem LC conforms to the DIN EN ISO 9917-2.

Keep the instructions for use for the duration of the application.

Order-No.	Article
7041215	15 g powder, shade Universal / 10 ml liquid dosage spoon, mixing pad